Bodenablassventil

Metall, DN 15 - 40

Tank Bottom Valve

Metal, DN 15 - 40

- Ⓓ ORIGINAL EINBAU- UND MONTAGEANLEITUNG
- ⒼB INSTALLATION, OPERATING AND MAINTENANCE INSTRUCTIONS

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Contents

1	General information	17
2	General safety information	17
2.1	Information for service and operating personnel	18
2.2	Warning notes	18
2.3	Symbols used	19
3	Definition of terms	19
4	Intended area of use	19
5	Technical data	20
6	Order data	20
7	Manufacturer's information	22
7.1	Transport	22
7.2	Delivery and performance	22
7.3	Storage	22
7.4	Tools required	22
8	Functional description	22
9	Construction	22
10	Installation and operation	23
10.1	Installing the valve	23
10.2	Operation	24
10.3	Shaft extension	24
11	Assembly / disassembly of spare parts	25
11.1	Valve disassembly (removing bonnet from body)	25
11.2	Removing the diaphragm	25
11.3	Mounting the diaphragm	25
11.3.1	General information	25
11.3.2	Mounting a concave diaphragm	26
11.3.3	Mounting a convex diaphragm	26
11.4	Bonnet mounting on the valve body	27
12	Commissioning	27
13	Inspection and servicing	28
14	Disassembly	28
15	Disposal	28
16	Returns	28
17	Information	29
18	Troubleshooting / Fault clearance	29
19	Sectional drawing and spare parts	30
20	EU declaration of conformity	31

1 General information

- Prerequisites to ensure that the GEMÜ valve functions correctly:
- x Correct transport and storage
 - x Installation and commissioning by trained personnel
 - x Operation according to these installation, operating and maintenance instructions
 - x Recommended maintenance
- Correct installation, operation, servicing and repair work ensure faultless valve operation.

	The descriptions and instructions apply to the standard versions. For special versions not described in these installation, operating and maintenance instructions the basic information contained herein applies in combination with any additional special documentation.
	All rights including copyright and industrial property rights are expressly reserved.

2 General safety information

- The safety information does not take into account:
- x Unexpected incidents and events, which may occur during installation, operation and servicing.
 - x Local safety regulations which must be adhered to by the operator and by any additional installation personnel.

2.1 Information for service and operating personnel

The installation, operating and maintenance instructions contain fundamental safety information that must be observed during commissioning, operation and servicing. Non-compliance with these instructions may cause:

- x Personal hazard due to electrical, mechanical and chemical effects.
- x Hazard to nearby equipment.
- x Failure of important functions.
- x Hazard to the environment due to the leakage of dangerous materials.

Prior to commissioning:

- Read the installation, operating and maintenance instructions.
- Provide adequate training for the installation and operating personnel.
- Ensure that the contents of the installation, operating and maintenance instructions have been fully understood by the responsible personnel.
- Define the areas of responsibility.

During operation:

- Keep the installation, operating and maintenance instructions available at the place of use.
- Observe the safety information.
- Use only in accordance with the specifications.
- Any servicing work and repairs not described in the installation, operating and maintenance instructions must not be performed without consulting the manufacturer first.

⚠ DANGER

Strictly observe the safety data sheets or the safety regulations that are valid for the media used.

In cases of uncertainty:

- x Consult the nearest GEMÜ sales office.

2.2 Warning notes

Wherever possible, warning notes are organised according to the following scheme:

⚠ SIGNAL WORD

Type and source of the danger

- Possible consequences of non-observance.
- Measures for avoiding danger.

Warning notes are always marked with a signal word and sometimes also with a symbol for the specific danger.

The following signal words and danger levels are used:

⚠ DANGER

Imminent danger!

- Non-observance will lead to death or severe injury.

⚠ WARNING

Potentially dangerous situation!

- Non-observance can cause death or severe injury.

⚠ CAUTION

Potentially dangerous situation!

- Non-observance can cause moderate to light injury.

CAUTION (WITHOUT SYMBOL)

Potentially dangerous situation!

- Non-observance can cause damage to property.

2.3 Symbols used

	Danger - hot surfaces!
	Danger - corrosive materials!
	Hand: indicates general information and recommendations.
	Bullet point: indicates the tasks to be performed.
	Arrow: indicates the response(s) to tasks.
	Enumeration sign

3 Definition of terms

Working medium

The medium that flows through the valve.

4 Intended area of use

- x The GEMÜ 643 tank bottom valve is welded into a tank bottom. It controls a flowing medium by manual operation.
- x **The valve may only be used providing the product technical criteria are complied with (see chapter 5 "Technical data").**
- x Do not paint the bolts and plastic parts of the valve!

⚠ WARNING

Use the valve only for the intended purpose!

- Otherwise the manufacturer liability and guarantee will be void.
- Use the valve only in accordance with the operating conditions specified in the contract documentation and in the installation, operating and maintenance instructions.
- The valve may only be used in potentially explosive zones confirmed in the declaration of conformity (ATEX).

5 Technical data

Working medium

Corrosive, inert, gaseous and liquid media which have no negative impact on the physical and chemical properties of the body and diaphragm material.

Maximum permissible pressure of working medium, applied upstream 10 bar

Temperatures

Media temperature -10 ... 100 °C

Sterilisation temperature ⁽¹⁾

EPDM (Code 13)	max. 150 °C ⁽²⁾ , max. 60 min per cycle
EPDM (Code 17)	max. 150 °C ⁽²⁾ , max. 180 min per cycle
PTFE/EPDM (Code 54)	max. 150 °C ⁽²⁾ , no time limit per cycle
PTFE/EPDM (Code 5M)	max. 150 °C ⁽²⁾ , no time limit per cycle


¹ The sterilisation temperature is valid for steam (saturated steam) or superheated water.

² If the sterilisation temperatures listed above are applied to the EPDM diaphragms for longer periods of time, the service life of the diaphragms will be reduced. In these cases, maintenance cycles must be adapted accordingly. This also applies to PTFE diaphragms exposed to high temperature fluctuations.

PTFE diaphragms can also be used as steam barriers; however, this will reduce their service life. The maintenance cycles must be adapted accordingly.

GEMÜ 555 and 505 globe valves are particularly suitable for use in the area of steam generation and distribution.

The following valve arrangement for interfaces between steam pipes and process pipes has proven itself over time: A globe valve for shutting off steam pipes and a diaphragm valve as an interface to the process pipes.

Ambient temperature 0 ... 60 °C

Nominal size [mm]	Bonnet Code	Operating pressure [bar]		Weight [kg]
		EPDM	PTFE	
15	2AT	0 - 10	0 - 10	3.0
20	2AT	0 - 10	0 - 10	3.0
25	2AT	0 - 10	0 - 10	3.0
32	3AT	0 - 10	0 - 10	6.0
40	3AT	0 - 10	0 - 10	6.0

6 Order data

Body configuration	Code
Tank valve body	B

Connection	Code
Butt weld spigots	
Spigots DIN	0
Spigots EN 10357 series B (formerly DIN 11850 series 1)	16
Spigot EN 10357 series A (formerly DIN 11850 series 2) / DIN 11866 series A	17
Spigots DIN 11850 series 3	18
Spigots SMS 3008	37
Spigot ASME BPE / DIN 11866 series C	59
Spigot ISO 1127 / EN 10357 series C / DIN 11866 series B	60
Spigots ANSI/ASME B36.19M Schedule 10s	63
Spigots ANSI/ASME B36.19M Schedule 5s	64
Spigots ANSI/ASME B36.19M Schedule 40s	65

Valve body material	Code
1.4435 (F316L), forged body	40
1.4435 (BN2), forged body Δ Fe<0.5 %	42

Diaphragm material	Code
EPDM	13
EPDM	17
EPDM	19
EPDM	36
PTFE/EPDM, two-piece	5M
Material complies with FDA requirements	

Control function	Code
Manually operated	0

Bonnet size	Code
Bonnet size 2AT (DN 15 - 25)	2AT
Bonnet size 3AT (DN 32 - 40)	3AT

Surface finish	Code
Code see page 21	

Internal surface finishes for forged bodies ¹

Readings for Process Contact Surfaces	Mechanically polished ²		Electropolished	
	Hygienic class DIN 11866	Code	Hygienic class DIN 11866	Code
Ra ≤ 0.80 µm	H3	1502	HE3	1503
Ra ≤ 0.60 µm	-	1507	-	1508
Ra ≤ 0.40 µm	H4	1536	HE4	1537
Ra ≤ 0.25 µm ³	H5	1527	HE5	1516

Readings for Process Contact Surfaces acc. to ASME BPE 2016 ⁴	Mechanically polished ²		Electropolished	
	ASME BPE Surface Designation	Code	ASME BPE Surface Designation	Code
Ra Max. = 0.76 µm (30 µinch)	SF3	SF3	-	-
Ra Max. = 0.64 µm (25 µinch)	SF2	SF2	SF6	SF6
Ra Max. = 0.51 µm (20 µinch)	SF1	SF1	SF5	SF5
Ra Max. = 0.38 µm (15 µinch)	-	-	SF4	SF4

¹ Surface finishes of customized valve bodies may be limited in special cases.

² Or any other finishing method that meets the Ra value (acc. to ASME BPE).

³ The smallest possible Ra finish for pipe connections with an internal pipe diameter < 6 mm is 0.38 µm.

⁴ When using these surfaces, the bodies are marked according to the specifications of ASME BPE.

The surfaces are only available for valve bodies which are made of materials (e.g. GEMÜ material codes 40, 41, 44) and use connections (e.g. GEMÜ connection codes 59, 80, 88) according to ASME BPE.

Ra acc. to DIN EN ISO 4288 and ASME B46.1

Order example	643	25	B	60	40	13	0	2AT	1503
Type	643								
Nominal size		25							
Body configuration (code)			B						
Connection (code)				60					
Valve body material (code)					40				
Diaphragm material (code)						13			
Control function (code)							0		
Bonnet size (code)								2AT	
Surface finish (code see page 4)									1503

7 Manufacturer's information

7.1 Transport

- Only transport the valve by suitable means. Do not drop. Handle carefully.
- Dispose of packing material according to relevant local or national disposal regulations / environmental protection laws.

7.2 Delivery and performance

- Check that all parts are present and check for any damage immediately upon receipt.
- The scope of delivery is apparent from the dispatch documents and the design from the order number.
- The performance of the valve is checked at the factory.

7.3 Storage

- Store the valve free from dust and moisture in its original packaging.
- Store the valve in "open" position.
- Avoid UV rays and direct sunlight.
- Maximum storage temperature: 40 °C.
- Solvents, chemicals, acids, fuels or similar fluids must not be stored in the same room as valves and their spare parts.


7.4 Tools required

- The tools required for installation and assembly are **not** included in the scope of delivery.
- Use appropriate, functional and safe tools.

8 Functional description

GEMÜ 643 is a 2/2-way metal tank bottom valve. It is manually operated with a side mounted gear and has an optical position indicator as standard. The stainless steel valve body is machined from a single block (no welds) and is designed for welding directly into a tank bottom. The distance piece and the gearbox housing are made of stainless steel. The bonnet can be rotated through 360° (see chapter 10.2 "Operation"). The tank bottom valve body is also available with a pneumatic or motorized actuator on request.

9 Construction

1 Valve body

2 Diaphragm

A Bonnet

10 Installation and operation

Prior to installation:

- Ensure that valve body and diaphragm material are appropriate and compatible to handle the working medium.
- **Check the suitability prior to the installation.**
See chapter 5 "Technical data".

10.1 Installing the valve

⚠ WARNING

The equipment is subject to pressure!

- Risk of severe injury or death!
- Only work on depressurized plant.

⚠ WARNING

Corrosive chemicals!

- Risk of caustic burns!
- Wear appropriate protective gear when installing.

⚠ CAUTION

Hot plant components!

- Risk of burns!
- Only work on plant that has cooled down.

⚠ CAUTION

Never use the valve as a step or an aid for climbing!

- This entails the risk of slipping-off or damaging the valve.

CAUTION

Do not exceed the maximum permissible pressure!

- Take precautionary measures to avoid possible pressure surges (water hammer).


- Installation work must only be performed by trained personnel.
- Use appropriate protective gear as specified in plant operator's guidelines.

Installation location:


⚠ CAUTION

- Do not apply external force to the valve.
- Choose the installation location so that the valve cannot be used as a foothold (climbing aid).
- Lay the pipeline so that the valve body is protected against transverse and bending forces, and also vibrations and tension.

- x Direction of the working medium: optional
- x Mounting position of the valve:

With shaft extension by user:

See also chapter 10.3 "Shaft extension".

Installation:

1. Ensure the suitability of the valve for each respective use. The valve must be appropriate for the piping system operating conditions (medium, medium concentration, temperature and pressure) and the prevailing ambient conditions. Check the technical data of the valve and the materials.
2. Shut off plant or plant component.
3. Secure against recommissioning.
4. Depressurize the plant or plant component.
5. Completely drain the plant (or plant component) and let it cool down until the temperature is below the media vaporization temperature and scalding can be ruled out.
6. Correctly decontaminate, rinse and ventilate the plant or plant component.

Installation - Butt weld spigots:

1. Adhere to good welding practices!
2. Disassemble the bonnet with the diaphragm before welding the valve body into the pipeline (see chapter 11.1).
3. Allow butt weld spigots to cool down.
4. Reassemble the valve body and the bonnet with diaphragm (see chapter 11.4).

Observe appropriate regulations for connections!


After the installation:

- Reactivate all safety and protective devices.

10.2 Operation


⚠ CAUTION	
	<p>Handwheel can become hot during operation!</p> <ul style="list-style-type: none"> ➤ Risk of burns! ● Ensure protective gloves are worn when operating handwheel.

Optical position indicator

Bonnet rotatable through 360°

Old version:


Loosen Allen setscrew (arrow) in the bonnet distance piece (do not remove). Rotate bonnet to desired position, fix position by tightening setscrew.

Allen key:

DN 15 - 25 / bonnet size 2	SW 2
DN 32 - 40 / bonnet size 3	SW 2.5

New version:

The upper and lower part of the bonnet are assembled with a ball bearing. Therefore the hexagon socket screw must not be removed to prevent the balls from falling out.

Loosen Allen setscrew (arrow) in the bonnet distance piece (**Do not remove**). Rotate bonnet to desired position, fix position by tightening setscrew.

Allen key:

DN 15 - 25 / bonnet size 2	SW 3
DN 32 - 40 / bonnet size 3	SW 3

10.3 Shaft extension




Important:


When the user installs a handwheel extension care should be taken that it has sufficient bearing points.

Torques for shaft extension:

DN 15 - 25 / bonnet size 2	11 Nm
DN 32 - 40 / bonnet size 3	14 Nm

11 Assembly / disassembly of spare parts

11.1 Valve disassembly (removing bonnet from body)

1. Move bonnet **A** to the open position.
2. Remove bonnet **A** from valve body **1**.
3. Move bonnet **A** to the closed position.

Important:

After disassembly, clean all parts of contamination (do not damage parts). Check parts for potential damage, replace if necessary (only use genuine parts from GEMÜ).

11.2 Removing the diaphragm

Important:

Before removing the diaphragm, please remove the bonnet, see "Valve disassembly (removing bonnet from body)".

1. Unscrew the diaphragm.
2. Clean all parts of the remains of product and contamination. Do not scratch or damage parts during cleaning!
3. Check all parts for potential damage.
4. Replace damaged parts (only use genuine parts from GEMÜ).

11.3 Mounting the diaphragm

11.3.1 General information

Important:

Mount the correct diaphragm that suits the valve (suitable for medium, medium concentration, temperature and pressure). The diaphragm is a wearing part. Check the technical condition and function of the diaphragm valve before commissioning and during the whole term of use. Carry out checks regularly and determine the check intervals in accordance with the conditions of use and / or the regulatory codes and provisions applicable for this application.

Important:


If the diaphragm is not screwed into the adapter far enough, the closing force is transmitted directly onto the diaphragm pin and not via the compressor. This will cause damage and early failure of the diaphragm and thus leakage of the valve. If the diaphragm is screwed in too far no perfect sealing at the valve seat will be achieved and the function of the valve is no longer ensured.

Important:


Incorrectly mounted diaphragm may cause valve leakage / emission of medium. In this case remove the diaphragm, check the complete valve and diaphragm and reassemble again proceeding as described above.

The compressor is loose.
Compressor and bonnet flange seen from below:

Place the compressor loosely on the bonnet spindle, fit the grooves **D** into the guides **C**. The compressor must be able to be moved freely between the guides!

11.3.2 Mounting a concave diaphragm


1. Move bonnet **A** to the closed position.
2. Place the compressor loosely on the bonnet spindle, fit the grooves **D** into the guides **C** (see chapter 11.3.1 "General information").
3. Check if the compressor fits closely in the guides.
4. Screw new diaphragm tightly into the compressor manually.
5. Check if the diaphragm boss fits closely in the recess of the compressor.

6. If it is difficult to screw it in, check the thread, replace damaged parts (only use genuine parts from GEMÜ).
7. When clear resistance is felt turn back the diaphragm anticlockwise until its bolt holes are in correct alignment with the bolt holes of the bonnet.


11.3.3 Mounting a convex diaphragm

1. Move bonnet **A** to the closed position.
2. Place the compressor loosely on the bonnet spindle, fit the grooves **D** into the guides **C** (see chapter 11.3.1 "General information").
3. Check if the compressor fits closely in the guides.
4. Invert the new diaphragm face manually; use a clean, padded mat with bigger nominal sizes.

5. Position the new backing diaphragm onto the compressor.
6. Position the diaphragm face onto the backing diaphragm.

- Screw diaphragm face tightly into the compressor manually. The diaphragm boss must fit closely in the recess of the compressor.

- If it is difficult to screw it in, check the thread, replace damaged parts.
- When clear resistance is felt turn back the diaphragm anticlockwise until its bolt holes are in correct alignment with the bolt holes of the bonnet.
- Press the diaphragm face tightly onto the backing diaphragm manually so that it returns to its original shape and fits closely on the backing diaphragm.

11.4 Bonnet mounting on the valve body

- Move bonnet **A** to the closed position.
- Open bonnet **A** approx. 50 %.
- Position bonnet **A** with the mounted diaphragm **2** on the valve body **1**, take care to align the compressor weir and valve body weir.
- Mount washers **19** and nuts **20**. Initially hand tighten only.
- Fully tighten the nuts **20** diagonally.
- Ensure that the diaphragm **2** is compressed evenly (approx. 10-15 %, visible by an even bulge to the outside).
- Check tightness of completely assembled valve.

Important:

Diaphragms set in the course of time. After valve installation and commissioning make sure to tighten the nuts **20** (see chapter 19 "Sectional drawing and spare parts").

12 Commissioning

⚠ WARNING

Corrosive chemicals!

- Risk of caustic burns!
- Check the tightness of the media connections prior to commissioning!
- Use only the appropriate protective gear when performing the tightness check.

⚠ CAUTION

Protect against leakage!

- Provide precautionary measures against exceeding the maximum permitted pressures caused by pressure surges (water hammer).

Prior to cleaning or commissioning the plant:

- Check the tightness and the function of the diaphragm valve (close and reopen the diaphragm valve).
- If the plant is new and after repairs rinse the piping system with a fully opened valve (to remove any harmful foreign matter).

Cleaning:

- x The plant operator is responsible for selecting the cleaning material and performing the procedure.

	<p>Important: Diaphragms set in the course of time. After valve installation and commissioning make sure to tighten the nuts 20 (see chapter 19 "Sectional drawing and spare parts").</p>
---	---

	<p>Important: Service and maintenance: Dependent on the operating conditions, regrease the threaded spindle. GEMÜ recommends the grease Boss-Fluorine Y 108/00 (99099484)</p>
--	--

13 Inspection and servicing

⚠ WARNING	
<p>The equipment is subject to pressure! ➤ Risk of severe injury or death! ● Only work on depressurized plant.</p>	

⚠ CAUTION	
	<p>Hot plant components! ➤ Risk of burns! ● Only work on plant that has cooled down.</p>

⚠ CAUTION	
<ul style="list-style-type: none"> ● Servicing and maintenance work may only be performed by trained personnel. ● GEMÜ shall assume no liability whatsoever for damages caused by improper handling or third-party actions. ● In case of doubt, contact GEMÜ before commissioning. 	

1. Use appropriate protective gear as specified in plant operator's guidelines.
2. Shut off plant or plant component.
3. Secure against recommissioning.
4. Depressurize the plant or plant component.

The operator must carry out regular visual examination of the valves dependent on the operating conditions and the potential danger in order to prevent leakage and damage. The valve also has to be disassembled in the corresponding intervals and checked for wear (see chapter 11 "Assembly / Disassembly of spare parts").

14 Disassembly

Disassembly is performed observing the same precautionary measures as for installation.

- Disassemble the valve (see chapter 11.1 "Valve disassembly (removing bonnet from body)").

15 Disposal

	<ul style="list-style-type: none"> ● All valve parts must be disposed of according to relevant local or national disposal regulations / environmental protection laws. ● Pay attention to adhered residual material and gas diffusion from penetrated media.
--	--

16 Returns

- Clean the valve.
- Request a goods return declaration form from GEMÜ.
- Returns must be made with a completed declaration of return.

If not completed, GEMÜ cannot process
 x credits or
 x repair work
 but will dispose of the goods at the operator's expense.

	<p>Note for returns: Legal regulations for the protection of the environment and personnel require that the completed and signed goods return declaration is included with the dispatch documents. Returned goods can be processed only when this declaration is completed.</p>
--	--

17 Information

	<p>Note on Directive 2014/34/EU (ATEX Directive): A supplement to Directive 2014/34/EU is included with the product if it was ordered according to ATEX.</p>
---	---

	<p>Note on staff training: Please contact us at the address on the last page for staff training information.</p>
---	---


Should there be any doubts or misunderstandings in the preceding text, the German version of this document is the authoritative document!

18 Troubleshooting / Fault clearance

Fault	Possible cause	Fault clearance
Medium escapes from leak detection hole*	Valve diaphragm faulty	Check valve diaphragm for damage, replace diaphragm if necessary
Valve doesn't open or doesn't open fully	Bonnet faulty	Replace bonnet
	Valve diaphragm incorrectly mounted	Remove bonnet, check diaphragm mounting, replace if necessary
Valve leaks downstream (doesn't close or doesn't close fully)	Operating pressure too high	Operate valve with operating pressure specified in data sheet
	Foreign matter between valve diaphragm and valve body weir	Remove bonnet, remove foreign matter, check valve diaphragm and valve body weir for damage and replace if necessary
	Valve body weir leaking or damaged	Check valve body weir for damage, if necessary replace valve body
	Valve diaphragm faulty	Check valve diaphragm for damage, replace diaphragm if necessary
Valve leaks between bonnet and valve body	Valve diaphragm incorrectly mounted	Remove bonnet, check diaphragm mounting, replace if necessary
	Bolting between valve body and bonnet loose	Retighten bolting between valve body and bonnet
	Valve diaphragm faulty	Check valve diaphragm for damage, replace diaphragm if necessary
	Valve body / bonnet damaged	Replace valve body / bonnet
Valve body connection to piping / tank leaks	Incorrect installation	Check installation of valve body in piping / tank
Valve body leaks	Valve body faulty or corroded	Check valve body for damage, replace valve body if necessary
Handwheel cannot be turned	Bonnet faulty	Replace bonnet
	Threaded spindle seized	Regrease the threaded spindle dependent on the operating conditions; replace bonnet if necessary (see chapter 11)

* see chapter 19 "Sectional drawing and spare parts"

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declare that the equipment listed below complies with the safety requirements of the Pressure Equipment Directive 2014/68/EU.

Description of the equipment - product type

Tank Bottom Valve
GEMÜ 643

Notified body: TÜV Rheinland Industrie Service GmbH
Number: 0035
Certificate no.: 01 202 926/Q-02 0036
Applied standards: AD 2000

Conformity assessment procedure:
Module H1

Note for equipment with a nominal size \leq DN 25:

The products are developed and produced according to GEMÜ process instructions and quality standards which comply with the requirements of ISO 9001 and of ISO 14001.

According to section 4, paragraph 3 of the Pressure Equipment Directive 2014/68/EU these products must not be identified by a CE-label.

Joachim Brien
Head of Technical Department

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Hereby we, **GEMÜ Gebr. Müller Apparatebau GmbH & Co. KG**